



LineSnuggler

By Nursing & Patient Care
Innovations LLC

Care Instructions:

Machine wash in hot water and
dry on high heat
every 3 days or when soiled

Use Instructions:

For single patient use. Discard after
1 month or patient discharge

Discard if any holes, rips, or stains
form. Do not sew or repair

Visit: www.LineSnugglers.com for use
instructions and to place an order
or request a sample

US FDA 513(g) reviewed:

Submission #C180053

Patent Pending - All rights reserved




Machine
Wash hot



Dry
High

MADE IN



U. S. A.



Lot 005